

Adam Equipment

ADAM THERMAL PRINTER

Software Revision V1.6AT & above



ENGLISH

Serial number of	the printer:	
Model of the pro	duct linked to this printer:	
Software revision printer:	number of the product linked to the	
Serial number of	the product linked to the printer:	
Date of Purchase:	:	
Name of the supp	plier and place:	

1.0 CONTENTS

(P.N. 3126611270 Revision B3, October 2016)

WARF	RANTY INFORMATION	Error! Bookmark not defined.6
9.	COMMAND FUNCTION LIST	11
8.	ATP BAUD RATE SETTING INSTRUCTIONS	10
7.	Serial Interface	9
6.	Front panel descriptions	
5.	Paper Roll Installation	
4.2	Printer connection	7
4.1	Packing List	
4.	Printer Installation and Operation	
3.	Technical Specifications	4
2.	SETTING UP	3
1.	INTRODUCTION	2



1. INTRODUCTION

- The Adam Thermal Printer is a general purpose thermal tally roll printer specifically designed to work with the Adam series of balances.
- The printer comes preset to work with the default RS-232 interface found on the Adam balances but can easily be changed to work on other equipment as well.
- The printer is easy to use with only an On/Off switch and a paper feed button.



2. SETTING UP

IMPORTANT

Clean the print head regularly to avoid poor print quality, only attempt to clean when the printer is off and the head has had time to cool. Clean with a soft cloth using diluted alcohol.

	Avoid extremes of temperature. Do not place in direct sunlight or near air conditioning vents.
	Make sure the scale is located on a strong table and free from vibration.
	Avoid unstable power sources. Do not use near large users of electricity such as welding equipment or large motors. Do not mix batteries and use only the factory approved power adapter supplied with the machine. Do not use batteries and the DC adapter at the same time.
	Keep free from vibration. Do not place near heavy or vibrating machinery.
	Avoid high humidity that might cause condensation. Keep away from direct contact with water. Do not spray or immerse the scales in water.
	Do not place near open windows, air-conditioning vents or fans that may cause a draft and unstable readings.
ACAUTION Electrostatic sensitive devices. To prevent equipment damage, use proper grounding techniques.	Avoid operating in high static areas, or weighing high static materials, such as plastics or powders, as this may damage electronics and affect measurements. Use grounding mats or bonding straps to reduce potential



3. TECHNICAL SPECIFICATIONS

Print method: Thermal line printing

• Paper width: 57.5 ± 0.5mm

• Effective print width: 48mm


• Print speed: 50mm/s or 14 lines/s

Print Head Performances:

Print columns: 384 dots/line or 8 dots/mm

Reliability print: 50Km

• Interface:



Serial interface: D-SUB 25 socket, optional baud rate 2400bps, 4800bps, 9600bps, 19200bps, 38400bps, 57600, 115200, default is at 4800bps, data structure no parity, 8 digital word size,1 digital stop bit, supports RTS/CTS and XON/XOFF handshake protocol.

Cash drawer: 6 wire RJ-11 socket, DC 12V/1A output signal to drive the cash drawer.

ΕN

• Characters:

ASCII character set: 12×24 dot, 1.25 (W) ×3.00 (H) mm

.

• Print paper:

Paper type: Thermal print paper

Paper width: 57.5±0.5mm,

Paper thickness: $0.06 \sim 0.08$ mm Roll outer dia. (max.): Φ 51mm

Recommend thermal print paper: ADAM Box of 10 – part number 3126011281 Box of 20 – part number 3126011282

• Printer power adaptor:

Adaptor input voltage: AC 110V/220V, 50~60Hz

Adaptor output voltage: DC 8.5V/3A

Printer weight: 0.8kg

Overall dimensions:

114 (W) × 184 (D) × 90 (H) mm

Working temperature: 0~45°C

Working relative humidity: 10∼80%



4. PRINTER INSTALLATION AND OPERATION


4.1 Packing List

Item	Part Number	Quantity
Adam Thermal Printer	2120011268	1
Power Adaptor Module	3124011269	1
Mains cable	Country dependant	1
RS 232 (25 pin to 9 pin)	3124011279	1
cable - female		
RS 232 (25 pin to 9 pin)	3124011273	1
cable - male		
Thermal paper	3126011263	1
Manual	3126611270	1

4.2 Printer connection


Please connect the cables as shown.

Make sure the printer and HOST machine are powered off and connect the data communication cable to the SERIAL PORT



The RJ-11 jack is for use when connecting the printer to a cash drawer.

Connect the Power adaptor to the printer and electricity supply (c)




Turn the printer on using the ON/OFF switch located at the rear of the printer

5. PAPER ROLL INSTALLATION

The thermal paper is fitted as below

- (1) Open the top cover
- (2) Place the paper roll into the space provided in the following way.



Correct direction

Wrong direction

- (3) Pull the leading edge of the paper outside of the printer and then close the top cover.
- (4) Tear off the excess paper using the cutter in the top cover

6. FRONT PANEL DESCRIPTIONS

FEED key: To feed the paper through

POWER light: Power light to indicate POWER is supplied to the printer

ERROR light: To indicate a paper jam or problem with printing


PAPER light: To indicate the printer is out of paper

7. SERIAL INTERFACE

Pin definition of serial interface:

Pins	Signal Description	Signal Source	Explanation
2	TXD	Printer	When you use XON/XOFF
			handshake protocol, printer
			sends control code XON/XOFF
			to the host machine.
3	RXD	Host machine	The printer receives data from
			the host machine
4	RTS	Printer	This signal indicates the printer
			status, low level indicates
			printer busy, can not receive
			data, the high level indicates
			printer ready, can receive data.
7	GND		Ground signal
20	DTR	Printer	Same with RTS signal(pin 4)

All Adam Equipment scales and balances use a 3 wire cable connection which includes TXD, RXD and GND. DTR and RTS are not required




Host machine serial interface DB-9



8. ATP BAUD RATE SETTING INSTRUCTIONS

Please contact Adam Equipment to obtain the application tool required

- 1. Connect the ATP to your PC using the cable provided
- 2. Power ON the ATP
- 3. Open the ATP baud rate setting application
- 4. Select the COM port that your PC is using
- 5. Select the current baud rate that your ATP is set to(default is 9600)
- 6. Select your required baud rate from the drop down list on the right hand side
- 7. Click on "Set Default Baud rate"
- 8. Click Close
- 9. Power your ATP OFF and then ON again





9. COMMAND FUNCTION LIST

No.	Comn	nands	Name
1	HT		Horizontal tab
2	LF		Print and line feed
3	CR		Print and carriage return
4	ESC	SO	Select Double Width
5	ESC	DC4	Disable Double Width
6	ESC	SP n	Set character right space amount
7	ESC	! n	Batch specify print mode
8	ESC	\$ nL nH	Specify absolute position
9	ESC	% n	Specify/cancel download character set
10	ESC	&	Define download characters
11	ESC	*	Specify bit image mode
12	ESC	- n	Specify/cancels underline mode
13	ESC	2	Set default line spacing
14	ESC	3 n	Set line feed amount
15	ESC	? n	Delete download characters
16	ESC	@	Initialize printer
17	ESC	D	Set horizontal tab position
18	ESC	En	Specify/cancel emphasized characters
19	ESC	Jn	Print and Paper Feed



WARRANTY STATEMENT

Adam Equipment offers Limited Warranty (Parts and Labour) for the components failed due to defects in materials or workmanship. Warranty starts from the date of delivery.

During the warranty period, should any repairs be necessary, the purchaser must inform its supplier or Adam Equipment Company. The company or its authorised Technician reserves the right to repair or replace the components at any of its workshops depending on the severity of the problems. However, any freight involved in sending the faulty units or parts to the service centre should be borne by the purchaser.

The warranty will cease to operate if the equipment is not returned in the original packaging and with correct documentation for a claim to be processed. All claims are at the sole discretion of Adam Equipment.


This warranty does not cover equipment where defects or poor performance is due to misuse, accidental damage, exposure to radioactive or corrosive materials, negligence, faulty installation, unauthorised modifications or attempted repair or failure to observe the requirements and recommendations as given in this User Manual. Additionally rechargeable batteries (where supplied) are not covered under warranty.

Repairs carried out under the warranty does not extend the warranty period.

Components removed during the warranty repairs become the company property.

The statutory right of the purchaser is not affected by this warranty. The terms of this warranty is governed by the UK law. For complete details on Warranty Information, see the terms and conditions of sale available on our web-site.

WEEE 2012/19/EU



This device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements. Disposal of batteries (if fitted) must conform to local laws and restrictions.

Cet appareil ne peut être éliminé avec les déchets ménagers. L'élimination de la batterie doit être effectuée conformément aux lois et restrictions locales.

Dieses Gerät nicht mit dem Hausmüll entsorgt.

Dispositivo no puede ser desechado junto con los residuos domésticos Dispositivo non può essere smaltito nei rifiuti domestici.

FCC / IC CLASS A DIGITAL DEVICE EMC VERIFICATION STATEMENT

NOTE: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules and Canadian ICES-003/NMB-003 regulation. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense..



Adam Equipment products have been tested with, and are always supplied with mains power adaptors which meet all legal requirements for the intended country or region of operation, including electrical safety, interference and energy efficiency. As we often update adaptor products to meet changing legislation it is not possible to refer to the exact model in this manual. Please contact us if you need specifications or safety information for your particular item. Do not attempt to connect or use an adaptor not supplied by us.



ADAM EQUIPMENT is an ISO 9001:2008 certified global company with more than 40 years' experience in the production and sale of electronic weighing equipment.

Adam products are predominantly designed for the Laboratory, Educational, Medical, retail and Industrial Segments. The product range can be described as follows:

- Analytical and Precision Balances
- Compact and Portable Balances
- High Capacity Balances
- Moisture analysers / balances
- Mechanical Scales
- Counting Scales
- Digital Weighing/Check-weighing Scales
- High performance Platform Scales
- Crane scales
- Medical Scales
- Retail Scales for Price computing

For a complete listing of all Adam products visit our website at

www.adamequipment.com

Adam Equipment Co. Ltd.	Adam Equipment Inc.	AE Adam GmbH
Maidstone Road, Kingston	1, Fox Hollow Rd.	Instenkamp 4,
Milton Keynes MK10 0BD	Oxford 06478	24242 Felde,
UK	USA	Deutschland.
Phone:+44 (0)1908 274545	Phone: +1 203 790 4774	T: +49- (0) 4340-40300-0
Fax: +44 (0)1908 641339	Fax: +1 203 792 3406	F: +49- (0) 4340-40300-20
e-mail: sales@adamequipment.co.uk	e-mail: sales@adamequipment.com	E: vertrieb@aeadam.de
Adam Equipment S.A. (Pty) Ltd.	Adam Equipment (S.E. ASIA) PTY	Adam Equipment (Wuhan) Co. Ltd.
7 Megawatt Road,	Ltd	A Building East Jianhua
Spartan EXT 22	2/71 Tacoma Circuit	Private Industrial Park
Kempton Park,	CANNING VALE 6155	Zhuanyang Avenue
Johannesburg,	Perth	Wuhan Economic & Technological
Republic of South Africa	Western Australia	Development Zone
Phone +27 (0)11 974 9745	Phone: +61 (0) 8 6461 6236	430056 Wuhan P.R.China
Fax: +27 (0)11 392 2587	Fax +61 (0) 8 9456 4462	Phone: + 86 (27) 59420391
e-mail: sales@adamequipment.co.za	E-mail:	Fax + 86 (27) 59420388
	sales@adamequipment.com.au	E-mail: info@adamequipment.com.cn

© Copyright by Adam Equipment Co. Ltd. All rights reserved. No part of this publication may be reprinted or translated in any form or by any means without the prior permission of Adam Equipment.

Adam Equipment reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

All information contained within this publication is to the best of our knowledge timely, complete and accurate when issued. However, we are not responsible for misinterpretations which may result from the reading of this material.

The latest version of this publication can be found on our Website www.adamequipment.co.uk